


Gewindedurchmesser	5,0mm
Kerndurchmesser	4,2mm
Kopfdurchmesser	6,8mm
Bohrer für das Gewinde	Ø 3,9mm
Schraubendreher	Ø 3,5mm

Winkelstabile Knochenschrauben 5

SSt	A
397 129 77 7981	16mm
397 129 77 7991	18mm
397 129 77 8001	20mm
397 129 77 8011	22mm
397 129 77 8021	24mm
397 129 77 8031	26mm
397 129 77 8041	28mm
397 129 77 8051	30mm
397 129 77 8061	32mm
397 129 77 8071	34mm
397 129 77 8081	36mm
397 129 77 8091	38mm
397 129 77 8101	40mm
397 129 77 8111	42mm
397 129 77 8121	44mm
397 129 77 8131	46mm
397 129 77 8141	48mm
397 129 77 8151	50mm
397 129 77 8161	55mm
397 129 77 8171	60mm
397 129 77 8181	65mm
397 129 77 8191	70mm
397 129 77 8201	75mm
397 129 78 7401	80mm
397 129 78 7411	85mm
397 129 78 7421	90mm
397 129 78 7431	95mm
397 129 78 7441	100mm
397 129 78 7451	105mm



Gewindedurchmesser	3,5mm
Kerndurchmesser	2,9mm
Kopfdurchmesser	6,8mm
Bohrer für das Gewinde	Ø 2,7mm
Schraubendreher	Ø 3,5mm

Winkelstabile Knochenschrauben 5/3,5

SSt	A
397 129 77 8231	16mm
397 129 77 8241	18mm
397 129 77 8251	20mm
397 129 77 8261	22mm
397 129 77 8271	24mm
397 129 77 8281	26mm
397 129 77 8291	28mm
397 129 77 8301	30mm
397 129 77 8311	32mm
397 129 77 8321	34mm
397 129 77 8331	36mm
397 129 77 8341	38mm
397 129 77 8351	40mm
397 129 77 8361	42mm
397 129 77 8371	44mm
397 129 77 8381	46mm
397 129 77 8391	48mm
397 129 77 8401	50mm
397 129 77 8411	55mm
397 129 77 8421	60mm
397 129 77 8431	65mm
397 129 77 8441	70mm
397 129 77 8451	75mm
397 129 78 7471	80mm
397 129 78 7481	85mm
397 129 78 7491	90mm
397 129 78 7501	95mm
397 129 78 7511	100mm
397 129 78 7521	105mm