



Gewindedurchmesser	4,5 mm
Kerndurchmesser	3,0 mm
Kopfdurchmesser	8,0 mm
Bohrer für das Gewinde	Ø 3,2 mm
Bohrer für Gleitöffnung	Ø 4,5 mm
Schraubendreher	Ø 3,5 mm

**Selbstschneidende Kortikalisschrauben HA 4,5**

SSt	A
<b>397 129 79 9421</b>	14 mm
<b>397 129 79 9431</b>	16 mm
<b>397 129 79 9441</b>	18 mm
<b>397 129 79 5521</b>	20 mm
<b>397 129 79 9451</b>	22 mm
<b>397 129 79 9461</b>	24 mm
<b>397 129 79 5541</b>	26 mm
<b>397 129 79 9471</b>	28 mm
<b>397 129 79 5551</b>	30 mm
<b>397 129 79 9481</b>	32 mm
<b>397 129 79 9491</b>	34 mm
<b>397 129 79 5571</b>	36 mm
<b>397 129 79 9501</b>	38 mm
<b>397 129 79 5581</b>	40 mm
<b>397 129 79 9511</b>	42 mm
<b>397 129 79 9521</b>	44 mm
<b>397 129 79 9531</b>	46 mm
<b>397 129 79 9541</b>	48 mm
<b>397 129 79 5601</b>	50 mm
<b>397 129 79 9551</b>	52 mm
<b>397 129 79 9561</b>	54 mm
<b>397 129 79 9571</b>	56 mm
<b>397 129 79 9581</b>	58 mm
<b>397 129 79 5621</b>	60 mm
<b>397 129 79 5631</b>	65 mm
<b>397 129 79 5641</b>	70 mm
<b>397 129 79 5651</b>	75 mm
<b>397 129 79 5661</b>	80 mm
<b>397 129 79 5671</b>	85 mm
<b>397 129 79 5681</b>	90 mm
<b>397 129 79 5691</b>	95 mm
<b>397 129 79 5701</b>	100 mm
<b>397 129 79 5711</b>	105 mm
<b>397 129 79 5721</b>	110 mm