



**ŠROUBY KOSTNÍ UZAMYKATELNÉ 5**  
 LOCKING BONE SCREWS 5


průměr závitu	thread diameter	5,0 mm
průměr jádra	core diameter	4,2 mm
průměr hlavy	head diameter	6,8 mm
vrták pro závit	drill bit for threaded hole	Ø 3,9 mm
šroubovák	screwdriver	Ø 3,5 mm

**Šroub kostní uzamykací 5**  
 Locking bone screw 5

SSt	Ti	A
397 129 77 7981	397 129 77 7984	16 mm
397 129 77 7991	397 129 77 7994	18 mm
397 129 77 8001	397 129 77 8004	20 mm
397 129 77 8011	397 129 77 8014	22 mm
397 129 77 8021	397 129 77 8024	24 mm
397 129 77 8031	397 129 77 8034	26 mm
397 129 77 8041	397 129 77 8044	28 mm
397 129 77 8051	397 129 77 8054	30 mm
397 129 77 8061	397 129 77 8064	32 mm
397 129 77 8071	397 129 77 8074	34 mm
397 129 77 8081	397 129 77 8084	36 mm
397 129 77 8091	397 129 77 8094	38 mm
397 129 77 8101	397 129 77 8104	40 mm
397 129 77 8111	397 129 77 8114	42 mm
397 129 77 8121	397 129 77 8124	44 mm
397 129 77 8131	397 129 77 8134	46 mm
397 129 77 8141	397 129 77 8144	48 mm
397 129 77 8151	397 129 77 8154	50 mm
397 129 77 8161	397 129 77 8164	55 mm
397 129 77 8171	397 129 77 8174	60 mm
397 129 77 8181	397 129 77 8184	65 mm
397 129 77 8191	397 129 77 8194	70 mm
397 129 77 8201	397 129 77 8204	75 mm
397 129 78 7401	397 129 78 7404	80 mm
397 129 78 7411	397 129 78 7414	85 mm
397 129 78 7421	397 129 78 7424	90 mm
397 129 78 7431	397 129 78 7434	95 mm
397 129 78 7441	397 129 78 7444	100 mm
397 129 78 7451	397 129 78 7454	105 mm



průměr závitu	thread diameter	3,5 mm
průměr jádra	core diameter	2,9 mm
průměr hlavy	head diameter	6,8 mm
vrták pro závit	drill bit for threaded hole	Ø 2,7 mm
šroubovák	screwdriver	Ø 3,5 mm

**Šroub kostní uzamykací 5/3,5**  
 Locking bone screw 5/3.5

SSt	Ti	A
397 129 77 8231	397 129 77 8234	16 mm
397 129 77 8241	397 129 77 8244	18 mm
397 129 77 8251	397 129 77 8254	20 mm
397 129 77 8261	397 129 77 8264	22 mm
397 129 77 8271	397 129 77 8274	24 mm
397 129 77 8281	397 129 77 8284	26 mm
397 129 77 8291	397 129 77 8294	28 mm
397 129 77 8301	397 129 77 8304	30 mm
397 129 77 8311	397 129 77 8314	32 mm
397 129 77 8321	397 129 77 8324	34 mm
397 129 77 8331	397 129 77 8334	36 mm
397 129 77 8341	397 129 77 8344	38 mm
397 129 77 8351	397 129 77 8354	40 mm
397 129 77 8361	397 129 77 8364	42 mm
397 129 77 8371	397 129 77 8374	44 mm
397 129 77 8381	397 129 77 8384	46 mm
397 129 77 8391	397 129 77 8394	48 mm
397 129 77 8401	397 129 77 8404	50 mm
397 129 77 8411	397 129 77 8414	55 mm
397 129 77 8421	397 129 77 8424	60 mm
397 129 77 8431	397 129 77 8434	65 mm
397 129 77 8441	397 129 77 8444	70 mm
397 129 77 8451	397 129 77 8454	75 mm
397 129 78 7471	397 129 78 7474	80 mm
397 129 78 7481	397 129 78 7484	85 mm
397 129 78 7491	397 129 78 7494	90 mm
397 129 78 7501	397 129 78 7504	95 mm
397 129 78 7511	397 129 78 7514	100 mm
397 129 78 7521	397 129 78 7524	105 mm